

Product Data

Product Name: Vaginal Speculum/Espéculo Vaginal

Size: S, M, L

Product Description

Vaginal speculums are composed of upon and down mobile parts along with its handle, and some with fixture on the handles. The fixtures can be fastener or screw with its nut

Intended Use

It can be used in checking and treatment gynecology disease in hospitals or Clinics

Specifications

Test Item		Technical Request
Physical Capacity 1	Back Outer Height h	Size should be measured up TKY/CE-TD-VS-03 3.1
	Appearance	Should be smooth ,no edges and corners on the cirque ,no burr, no dust and no watermark
	Close neatly and no baffling	Vaginal Speculum should operate with single hand, close neatly and no baffling
	Axes and axes hole	Operate neatly between axes and axes holes, no casting when used
	Expanding neatly	Expand neatly and no baffling
	No parts prevent eyeshot	No parts prevent the eyeshot
	Flexibility	Vaginal Speculum expands the last position , the flexibility value $\leq 30\text{mm}$ under the 20N rated loading
Physical Capacity 2	Intensity	Vaginal Speculum expands the last position, close to the upper part 20mm,any part should not rupture under 30N quiet loading
	l	Measured up TKY/CE-TD-VS-03 3.1
	w	Measured upTKY/CE-TD-VS-03 3.1
	b	Measured up TKY/CE-TD-VS-03 3.1
Biology Capability	Top Difference between Upper parts and lower parts	Top Difference between Upper parts and lower parts $\leq 2\text{mm}$
	Sterile	Should be sterile
	Cell toxicity	The reaction is no more than 1 grade.
	Sensitization	The skin reaction is no more than I grade.
	Skin irritatin	No skin reaction, P II is 0~ 0.4.

Product Picture

Vaginal Speculum Lateral Screw Type



L(white nut) M(red nut) S(yellow nut)

Vaginal Speculum Middle Screw Type



L M/L M S

Vaginal Speculum Fastener Type



L(green fastener) M(red fastener) S(blue fastener) XS(white fastener)

Vaginal Speculum Push Type



Large Medium Small

Vaginal Speculum Fast Blocking Type



Large Medium Small

Vaginal Speculum Spanish Type



Large Medium Small